



Kingfisher

Solution & Blanket

Warming Cabinets



Designed specifically for use in hospitals, the Kingfisher range of Solution and Blanket Warming Cabinets provide an accurate and cost-effective means for storing and warming a wide variety of clinical items such as:

- Infusion fluids
- Filled syringes
- Lotions
- Irrigation fluids
- Ointments
- Blankets

Kingfisher Solution Warming Cabinets

When designing the Kingfisher range, special consideration was given to the correct storage of bags and bottles. The special rack for IV bags prevents overloading of the cabinets and allows adequate air circulation around each bag, thus maximising the cabinet's performance. The control panel has been designed to help prevent accidental or unauthorised tampering with the temperature control or ON/OFF switch. Temperatures are PID controlled and a digital temperature display constantly monitors chamber conditions.

The Kingfisher range of Solution Warming Cabinets has been designed to accept a certain number of 1 litre bottles and/or bags.

The table below gives an indication of the possible storage capacities for each model. Please note that with every model except LK/SWC/01/1, a mixture of bottles and storage racks can be accommodated.

Kingfisher Blanket Warming Cabinets are suitable for storing solutions as well as blankets and are supplied with 2 or 4 shelves, depending on the model. These cabinets are not fan circulated, ensuring maximum storage and safety.

NB – Blankets must never be stored in a Solution Warming Cabinet.

Features in the range include:

- Precise electronic temperature control
- Tamper-proof controls
- Over-temperature protection
- Radiused stainless steel interior
- Fan circulated (Solution Warming Cabinets only)
- Specially designed racks for I.V. bags (optional)

Specifications

Solution Warming Cabinets

Cat. No.	Nominal capacity (Litres)	Internal Dimensions H x W x D mm	External Dimensions H x W x D mm	No. of shelves	Power rating (Watts)	Air circulation	Temperature range °C	Storage capacity (1 Litre bottles)	Rack capacity (12 x IV bags per rack)
LK/SWC/01/1	72	350 x 520 x 400	570 x 620 x 630	1	250	Fan circulation	5 above ambient to 80	12	1 (12)
LK/SWC/02/1	145	700 x 520 x 400	920 x 620 x 630	2	350			24	2 (24)
LK/SWC/03/1	176	850 x 520 x 400	1070 x 620 x 630	3	550			36	2 (24)*
LK/SWC/04/1	228	1100 x 520 x 400	1320 x 620 x 630	4	750			48	3 (36)*

*These models are designed with spare storage space at the top for smaller bottles or other items

Blanket Warming Cabinets

Cat. No.	Nominal capacity, (Litres)	Internal Dimensions H x W x D mm	External Dimensions H x W x D mm	No. of shelves	Power rating (Watts)	Air circulation	Temperature range °C
LK/BWC/02/1	145	700 x 520 x 400	1100 x 620 x 630	2	350	Natural convection	5 above ambient to 50
LK/BWC/04/1	228	1100 x 520 x 400	1320 x 620 x 630	4	750		

Options and Accessories

- LK/SWC/10/0** Storage rack to accommodate 12 x 1 litre bags.
LK/INC/11/0 Spare shelf
LK/BWC/12/0 Spare shelf

LTE Scientific Ltd.

Greenbridge Lane, Greenfield
Oldham, OL3 7EN, United Kingdom.

Tel: +44 (0) 1457 876221 • Fax: +44 (0) 1457 870131
Email: info@lte-scientific.com • Website: www.lte-scientific.com

